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OPERATING INSTRUCTIONS

PRODUKTBESCHREIBUNG

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OPERATING INSTRUCTIONS

PRODUCT DESCRIPTION

The use of far-infrared heat energy is a proven method of complementing beauty and wellness routines and naturally supporting the appearance of the skin.

Far infrared radiation works evenly without overheating the skin. It generates a gentle, penetrating warmth that envelops the body, promoting relaxation and contributing to increased well-being.

The pleasant feeling of warmth contributes to relaxation and promotes a natural radiance – for more glow and elasticity.

The system uses modern, evenly distributed far-infrared heating elements and an intelligent temperature control system that keeps the heat at a constant, comfortable level. This makes the application a soothing self-care moment – for visibly well-groomed skin and an all-round relaxed body feeling.

INTENDED USE

This device is intended to promote relaxation and create a pleasant feeling of warmth as part of wellness and beauty treatments.

It can be used as part of a personal relaxation or self-care routine to support general well-being and improve the feel of the skin.

The device is not intended for medical purposes or for the diagnosis, alleviation, or treatment of diseases.

SAFETY INSTRUCTIONS AND WARNINGS

It is important that you read all warnings and precautions in this manual. They are for your safety, to prevent injury, and to avoid damage to the device.

- Only use the device as described in this manual.
- **Do not use if you are pregnant**, wearing a pacemaker, have severe cardiovascular disease, acute skin injuries, or inflamed skin.
- Not suitable for children or persons with impaired perception, except under supervision.
- Do not use in damp environments (bathroom, shower).
- Do not allow liquids to enter the device.
- Always insert the power plug completely into a properly installed outlet.
- Do not cover the device during operation (except with the textile covers provided).
- If the cable, plug, or device is damaged, unplug it immediately.
- Keep away from flammable materials.
- Only clean the device when it is switched off and has cooled down.

PRECAUTIONS

- Before each use, connect the 4-pin plug of the arm/leg cuffs correctly to the front of the control unit.
- Only then switch on the device at the main switch (rear).
- After each session, first switch off the device using the **POWER button**, then press the main switch on the rear.
- When the arm/leg cuffs are empty, the temperature may briefly exceed the set value. As soon as the cuff is in contact with the body, the temperature will regulate itself to the set value with only a slight deviation. This is normal.
- Clean all textiles thoroughly after each use to avoid residues and extend the service life of the device.
- If the device is to be used in conjunction with other cosmetic devices or slimming aids, it is essential to consult a doctor beforehand.

SCOPE OF DELIVERY

- 1× main device
- 1× infrared leg cuffs
- 2× infrared arm cuffs
- 1× power cord
- 1× operating instructions

TECHNICAL DATA

- Model: PH-2A
- Supply voltage: AC 220–240 V~ 50 Hz
- Output voltage: DC 36 V
- Power consumption: max. 300 W
- Temperature range: 30–60 °C (standard 42 °C)
- Timer: 5–95 min. (standard 40 min.)
- Channels: 4, individually controllable
- Protection class: II (double insulation)
- Overload protection: yes
- Error display: yes
- Material: 100% polyester (PES) bonded with polyvinyl chloride fibers (CLF)

OPERATING INSTRUCTIONS

1. Connect

- Connect the textile covers (trouser and arm cuffs) to the corresponding sockets via the 4-pin plugs (note the channels).
- Connect the power cord and turn on the main switch (on the back).

2. Start

- Press the „POWER“ button.
- Select the desired channel with „ZONE“.
- Set the temperature (30–60 °C) using „Temperature +/-“.
- Set the time (5–95 min) using „Time +/-“.
- Start or stop heating using „ON/OFF“.

3. Switching off

- Press the „POWER“ button.
- Switch off the main switch and pull out the mains plug.

SYSTEM STRUCTURE & FUNCTION

- **Control unit:** regulates temperature and time for each zone.
- **Infrared leg cuffs:** gentle warmth for the waist, hips, and legs.
- **Infrared arm cuffs:** warmth for the arms.
- **Display & controls:** Temperature, time, channel status display; easy to set.
- **Safety features:** automatic overheating protection and error display in case of malfunctions.

MAINTENANCE

- Clean textile covers after each use with a soft, slightly damp cloth.
- Do not use aggressive cleaners or solvents.
- Wipe the control unit with a dry cloth only.

STORAGE & TRANSPORT

- Store in a cool, dry place away from direct sunlight.
- Roll up the cable loosely; do not kink or crush it.
- Transport only in the original packaging or equivalent protection.

TROUBLESHOOTING

- E0 – Overload → Reduce load or lower temperature.
- E1–E4 – Sensor circuit open (channels 1–4) → Check connection and fabric cover.
- E5–E8 – Sensor short circuit (channel 1–4) → Switch off the device and contact service.

SYMBOLS & EXPLANATIONS

	Device complies with EU regulations
	WEEE – Do not dispose of in household waste
	Read the operating instructions
	Manufacturer
	Manufactured in: China
	Series/batch number
	Packaging
	France - Disposal instructions for electrical appliances
	France - Disposal instructions for packaging
	Do not fold or crease.
	No sharp objects
	Not suitable for children under 3 years of age
	Do not wash
	Do not bleach
	No dry cleaning

MANUFACTURER/IMPORTER INFORMATION

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DISPOSAL

GERMANY:

Electronic products must not be disposed of in unsorted municipal waste. Instead, use the available return and collection systems. In this way, we can all contribute to the reuse, recycling, and recovery of electronic products.

Incorrect disposal of electronic products not only wastes resources (most electrical appliances contain valuable raw materials such as gold, copper, silver, etc.), but also endangers the environment due to the sometimes hazardous substances they contain.

ALTERNATIVE: We are happy to dispose of your old Reboots free of charge. Simply send them to:

WELLCOSAN GmbH
Wellcosan GmbH
Bellinger Tor 16a
36396 Steinau a. d. Str.

Please include a note requesting disposal.

FRANCE:

Electronic products must not be disposed of in unsorted municipal waste. Instead, use the available return and collection systems. In this way, we can all contribute to the reuse, recycling, and recovery of electronic products.

Incorrect disposal of electronic products not only wastes resources (most electrical appliances contain valuable raw materials such as gold, copper, silver, etc.), but also endangers the environment due to the sometimes hazardous substances they contain.

SLOVAKIA:

Electronic products must not be disposed of in unsorted municipal waste. Instead, use the available return and collection systems. In this way, we can all contribute to the reuse, recycling, and recovery of electronic products.

Improper disposal of electronic products not only wastes resources (most electrical appliances contain valuable raw materials such as gold, copper, silver, etc.), but also endangers the environment due to the sometimes hazardous substances they contain.

IMPORTANT INFORMATION ON ELECTROMAGNETIC COMPATIBILITY (EMC)

1. This product requires special precautions regarding EMC and must be installed and operated in accordance with the EMC information provided. This device may be affected by portable and mobile RF communications equipment.
2. Do not use a mobile phone or other devices that emit electromagnetic fields near the device. This may cause the device to malfunction.
3. Caution! This device has been thoroughly tested and inspected to ensure proper performance and function!
4. Caution: This device should not be used next to or stacked with other devices. If use next to or stacked with other devices is necessary, this device should be monitored to verify normal operation in the configuration in which it is used.

WARRANTY STATEMENT

We provide a 24-month warranty on the control unit. In the event of a malfunction, it may be necessary to send the device in for inspection. Please ensure that the shipping box is adequately padded to prevent any damage during transport. Unfortunately, no warranty claims can be made for defects that are attributable to improper shipping. Warranty conditions: Insofar as there is a legal obligation to provide a warranty, either a replacement will be provided free of charge or the manufacturer will be given the opportunity to repair the device, excluding the right to rescission or reduction. If, despite several attempts, the repair is still defective or a replacement device is also defective for reasons for which the manufacturer is responsible, the customer is entitled to rescission or reduction.

Warranty exclusion: The warranty does not cover defects caused by violent damage, improper operation, external force, or modifications and repairs carried out by third parties, such as defects caused by incorrectly dimensioned or short-circuited fuses, or defects attributable to normal wear and tear.

Warranty requirement: The warranty claim is only valid in conjunction with the original proof of purchase (invoice). Therefore, please keep the proof of purchase in a safe place. If you have any problems or questions when using the system, please do not hesitate to contact us at +49 (0) 157 359 900 84 or service@venenengel.de.

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Wellcosan GmbH c/o Caya Postbox 794627